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Referinte: Matopat.ro

Scutece-adulti.ro

	dry, black necrosis	yellow necrosis	red granulating tissue	pink epithelising wound	yellow infected necrolysis					
HEALING PHASE										
COLOURFUL WOUND CLASSIFICATION MATRIX										
CHALLENGE	dissolving the dry necrotic tissue initiate the healing process	debridement of necrotic tissue	promoting the growth of granulation tissue	stimulation of epithelisation process protection of migrating epidermal	protection against infection inhibition of infection					
RECOMMENDED PROCEEDINGS	surgical excision and / or stimulating the process of autolysis maintaining a moist environment	maintaining a moist environment absorbing excess exudate and necrotic material	maintaining a moist environment absorbing excess exudate protecting against mechanical injury	maintaining a moist environment stimulating epidermal growth protecting against injury	inhibiting the multiplication of bacteria					
CHOICE OF THE PROPER DRESSING										
EXUDATE										
SUITABLE DRESSING	Medisorb G	Medisorb A Medisorb P Medisorb P PLUS	Medisorb A Medisorb P Medisorb H Medisorb P PLUS	Medisorb H	Medisorb A Medisorb P Medisorb P PLUS	Medisorb A Medisorb P Medisorb H Medisorb P PLUS	Medisorb H	Medisorb F	Medisorb SILVER PAD	Medisorb SILVER
Level of exudate:										